


# Application Note

**Comparison between  
Precipitation & Size Exclusion  
Chromatography (SEC)  
vs. Membrane Affinity  
for exosome isolation**



# Exosome isolation from cell culture conditioned medium using Exo-spin™ Midi Columns and Membrane Affinity

Data Courtesy of Dr. Triantafillos Liloglou and Dr. Amelia Acha-Sagredo  
University of Liverpool, United Kingdom (Issued March, 18<sup>th</sup> 2019)

- **Experiment**  
Isolate exosomes using Exo-spin™
- **Exosomes origin**  
Conditioned cell culture medium
- **Initial sample volume**  
10 ml per Exo-spin™ Midi Column
- **Elution sample volume**  
500 µl per fraction in PBS

---

## Summary

In an effort to reduce lung cancer mortality rates through understanding mechanisms and improved early diagnosis, much research focus has recently been oriented towards extracellular vesicles (EVs). This is both in the context of EV biogenesis and also biomarker potential.

The research reported here compared two exosome isolation methods applied to conditioned cell culture medium from lung cancer cells. Results from samples pre-treated with Exo-spin™ Buffer and purified using Exo-spin™ Midi Columns were compared to the purification efficacy of a membrane affinity spin column method.

The Exo-spin™ purification system has been identified as superior to the membrane affinity technique. This was described through two major benefits: 1) exosomes do not aggregate; and 2) exosomes are ready to use without further buffer exchange, even in functional downstream applications. Analyses such as particle number and vesicle size have been performed.

## Methods


### Exo-spin™ Midi Columns

1. Remove cells and cellular debris
2. Precipitate exosomes containing pellet using Exo-spin™ Buffer (optional)
3. Purify exosomes using Exo-spin™ Midi SEC Columns

### Membrane affinity



1. Remove cells and cellular debris
2. Mix samples with binding buffer and bind to column
3. Remove unbound particles by washing
4. Elute EVs in elution buffer

## Results



**(A) Comparison of total yields from 10 ml conditioned cell culture medium by EV isolation method.** Total particle yield obtained with Exo-spin™ Buffer + Exo-spin™ Midi Columns is greater than the yield obtained using the membrane affinity method. Notably, Exo-spin™ Buffer is much more efficient than generic PEG 6000 when used prior to Exo-spin™ Midi Column exosome isolation.

**(B) NTA of EV fractions.** EV samples were isolated from 10 ml conditioned medium. **(1)** The overwhelming majority of particles isolated with the Exo-spin™ technology were in the 100-150 nm range, following a narrow size distribution expected for highly pure, high quality exosomes. **(2)** Vesicles isolated using the membrane affinity method, showed a broad range of particles with the majority in the 200-400 nm region. Blue numbers indicate peak particle sizes.



© 2013-2024 Cell Guidance Systems. All rights reserved. The trademarks mentioned herein are the property of Cell Guidance Systems or their respective owners.

Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*

#### Growth Factors

- Conventional (unformulated)
- PODS® - Sustained release

#### Exosomes

- Exo-spin™ - Purification
- ExoLISA™ - ELISA-like detection
- Instant Exosomes™ - purified and characterized
- NTA Service
- Freeze drying service

#### PeptiGel®

- Tunable self-assembling peptide hydrogels

#### Other products and services

- Small Molecules
- Softwell™ - 2D hydrogel (Europe only)
- Orangu™ - Cell counting reagent
- LipoQ™ - Lipid quantification assay
- Primary Hepatocytes

#### Cytogenetics

- Karyotype Analysis
- Array Hybridization

Scan for Exo-spin product page



General [info@cellgs.com](mailto:info@cellgs.com)  
Technical Enquiries [tech@cellgs.com](mailto:tech@cellgs.com)  
Orders and Quotes [order@cellgs.com](mailto:order@cellgs.com)

[www.cellgs.com](http://www.cellgs.com)

**EUROPE**  
**Cell Guidance Systems Ltd**  
Maia Building  
Babraham Bioscience Campus  
Cambridge  
CB22 3AT  
United Kingdom  
T +44 (0) 1223 967316  
F +44 (0) 1223 750186

**USA**  
**Cell Guidance Systems LLC**  
Helix Center  
1100 Corporate Square Drive  
St. Louis  
MO 63132  
USA  
T 760 450 4304  
F 314 485 5424