

Application Note

**Exosome marker detection:
an exquisitely sensitive
Europium Time Resolved
Immunofluorescence
assay**



Exosome detection from cell culture media using ExoLISA™ assay and BMG LABTECH TRF microplate readers

Data Courtesy of Laur-Alexandru Botos and Manoja Rasamanikkam
Cell Guidance Systems and [BMG LABTECH](#), United Kingdom (Issued August, 5th 2019)

- **Experiment**
Detect exosome markers using ExoLISA™ assay combined with BMG LABTECH plate readers
- **Exosomes origin**
25 µg exosomes isolated from conditioned cell culture media (prostate and colorectal cell lines) using Exo-spin™ exosome purification kit
- **Concentration range**
Serial dilutions in PBS from 100 µg/ml to 0.8 µg/ml, obtaining a total of 8 different exosome concentrations
- **Volume used per well**
100 µl of each dilution

Summary


The rapidly expanding EV research field benefits tremendously from technologies which enable the detection of exosomes and their content. Tetraspanin proteins such as CD9, CD63, and CD81 with exposed domains are particularly enriched in the exosome membrane and are often used as exosome biomarkers. ExoLISA™ - the Time Resolved Immunofluorescence exosome detection assay - has been developed specifically for analyzing the levels of these tetraspanin proteins on the surface of exosomes. ExoLISA™ assay enables researchers to directly measure tetraspanins present on exosomes from sources such as cell culture media, urine, or blood with low background signals and high sensitivity.

The assay has been successfully evaluated on several BMG LABTECH plate readers, including FLUOstar Omega, CLARIOstar Plus and the PHERAstar FSX. Accurate and reliable TRF signals were obtained for different exosome markers, with all instruments meeting the sensitivity requirements, ensuring optimal assay performance. Consistently linear signals were generated over a wide range of exosome concentrations, demonstrating the robustness of the assay. The ExoLISA™ assay can also be used with unpurified samples.

- **Detects different expression levels** of exosome markers
- **Simple** assay offering a high degree of reproducibility

- **Highly sensitive** and specific: only antigens displayed in multiple copies are detectable

Methods



Step 1: Biotinylated antibody is bound to streptavidin coated assay plates. The kit provides a 96-well microplate with an 8-well strip format.



Step 2: Biological samples containing exosomes are added to the plate. Exosomes bearing the appropriate antigen (CD9, CD63, or CD81) and any free antigen are captured by the immobilized antibody.

Step 3: Europium-labeled antibody is added and binds specifically to exosome antigens. Free antigen is not bound by the europium-labeled antibody (since the epitope is already occupied by the immobilized antibody) and is therefore not detected. Samples are then read on a Time-Resolved Fluorescence (TRF) plate reader such as the PHERAstar FSX.

Table 1. Instrument settings

Optic settings	Time resolved fluorescence (TRF)	
	Filters	Excitation: 337 nm Emission: 615 nm
General settings	Number of flashes	200 with flash lamp
	Settling time	0.3 s
	Integration start	200 μ s
	Integration time	400 μ s

Results



(A) Marker detection on exosomes from prostate and colorectal cancer cells, respectively. Each datapoint represents an average value (N=3). The error bars for each datapoint are very small, indicating extremely little variance in signal output. Moreover, the R² values for each serial dilution demonstrate excellent linearity. This is consistent with the data obtained for total protein concentration.

(B) Signal to blank ratio of CD81/CD63 detection on exosomes from prostate and colorectal cancer cells, respectively. Each datapoint represents an average value (N=3). The data presented further demonstrates the consistency of the TRF signal generated by the assay over a wide range of concentrations.

References

- Marina Colombo *et al.* Annu. Rev. Cell Dev. Biol. (2014) 30: 255-89. DOI: 10.1146/annurev-cellbio-101512-122326
- Diedrik Duijvesz *et al.* Int. J. Cancer (2015) 137: 2869-78. DOI: 10.1002/ijc.29664
- www.bmglabtech.com

For more information on our exosome range, please visit our website www.cellgs.com.



Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*

Growth Factors

- Recombinant
- Sustained Release

Exosomes

- Purification
- Detection
- NTA Service

Matrix Proteins

Small Molecules

Cytogenetics Analysis

Cell Counting Reagent



General info@cellgs.com
 Technical Enquiries tech@cellgs.com
 Quotes quotes@cellgs.com
 Orders order@cellgs.com

www.cellgs.com

EUROPE

Cell Guidance Systems Ltd
 Maia Building
 Babraham Bioscience Campus
 Cambridge
 CB22 3AT
 United Kingdom
 T +44 (0) 1223 967316
 F +44 (0) 1223 750186

USA

Cell Guidance Systems LLC
 Helix Center
 1100 Corporate Square Drive
 St. Louis
 MO 63132
 USA
 T 760 450 4304
 F 314 485 5424