


# Application Note

**Modulating protein  
release rates from  
PODS<sup>®</sup> crystals**



# Pre-incubating PODS<sup>®</sup> crystals to modulate release of GM-CSF

## Introduction to PODS<sup>®</sup>

### The challenge with soluble growth factors

Many proteins, especially growth factors and cytokines, when used as a reagent, degrade quickly, rapidly losing their bioactivity. This fragility hampers research and significantly limits the therapeutic potential of proteins.

### Protein Micro-depots


Development of a technology that can continuously replenish active protein from a local, microscopic store has been a significant challenge, but one that could transform the fields of cell culture and medicine by allowing greater control over the growth of cells.

### Introducing PODS<sup>®</sup>

PODS<sup>®</sup> technology has made the goal of a micro-depot for proteins a reality. PODS<sup>®</sup> is a sustained release system which continuously replenishes proteins from millions of local microscopic stores which can be placed next to (or at a distance from) cells, either randomly or in precise locations. Just like cells, these micro-depots release a steady stream of bioactive protein. This protein can be limited to local surroundings or dispersed more widely, or made to form a gradient.

### How does it work?

At the heart of PODS<sup>®</sup> is an extraordinary polyhedrin protein. This specific polyhedrin protein has the unique ability to encase cargo proteins within perfect, transparent, cubic, micro-sized crystals, much smaller than the cells. These protein crystals form admixtures of the polyhedrin and cargo proteins which slowly degrade releasing the biologically active cargo protein.



### How can PODS<sup>®</sup> help my research?

PODS<sup>®</sup> are tough and will withstand physical and chemical stress, so you can handle them with ease. PODS<sup>®</sup> can be made to release intact cargo protein over days, weeks or even months. Using PODS<sup>®</sup> you can readily create a steady-state protein environment in microscopic detail wherever you want, tailored exactly to your requirements. This is the power of PODS<sup>®</sup>. PODS<sup>®</sup> proteins are now available for many growth factors and cytokines and are already being used in many leading world-class research labs. PODS<sup>®</sup> protein applications include:

- Micropatterning
- Physiological, stable gradient formation
- Bioinks for 3D printing
- Microcarriers
- Functionalizing scaffolds
- Microfluidics (lab on a chip)
- Improved and simplified stem cell culture
- Therapeutic protein delivery


## Methods

**Pre-incubation:** PODS<sup>®</sup> GM-CSF crystals and PODS<sup>®</sup> Empty crystals ( $5 \times 10^5$ ) were spotted into wells of a 96-well plate and dried on. Subsequently, RPMI + 10% FBS was added to each well and incubated at 37°C.

**Culture method:** PODS<sup>®</sup> GM-CSF crystals and PODS<sup>®</sup> Empty crystals were spotted into wells as described above. TF-1 cells, which are dependent on GM-CSF, were then seeded and cultured for 5 days. **NOTE:** a single application of PODS<sup>®</sup> crystals was used during the culture period without any media change.


**Cell counting:** TF-1 cell number was assessed using the cell counting solution Orangu<sup>™</sup> (Cat OR01) according to the guidelines. Briefly, 10  $\mu$ l of Orangu<sup>™</sup> solution per 100  $\mu$ l of cell culture medium was added to each well. The plate was then incubated at 37°C for 2 hours, and subsequently centrifuged at 3000  $xg$  for 20 minutes to prevent carry-over of PODS<sup>®</sup> crystals. The supernatant was transferred into wells of a fresh plate and the absorbance measured at 450 nm using a microplate reader.

## Results



**Release of GM-CSF over 24 hours from PODS<sup>®</sup> GM-CSF crystals, quantified by ELISA.** PODS<sup>®</sup> GM-CSF crystals ( $5 \times 10^5$ ) were spotted onto 96-well plates and dried on. Subsequently, RPMI + 10% FBS was added to each well and incubated at 37°C. Medium was removed at indicated time points. GM-CSF was quantified by ELISA. Error bars represent 3 technical repeats.

**Proliferation of TF-1 cells in the presence of PODS<sup>®</sup> GM-CSF with or without pre-incubation.** PODS<sup>®</sup> GM-CSF crystals or PODS<sup>®</sup> Empty crystals were spotted onto a 96-well plate, after which RPMI + 10% FBS was added and incubated for 24 h. Subsequently,  $2 \times 10^3$  TF-1 cells in RPMI + 10% FBS were directly seeded on top and incubated for a further 5 days (blue bars). Cell number was assessed using a colorimetric assay, and proliferation was plotted relative to unsupplemented TF-1 cells. Error bars represent 8 technical repeats.



## Conclusions

- Serum-containing cell culture medium can activate the release of cargo protein from PODS<sup>®</sup> crystals.
- Pre-incubating PODS<sup>®</sup> crystals provides a starting amount of cargo protein in culture medium, beneficial if an initial supply of protein is critical.
- A single application of PODS<sup>®</sup> crystals is effective, significantly reducing both hands-on time and cost of materials.

© 2013-2024 Cell Guidance Systems. All rights reserved. The trademarks mentioned herein are the property of Cell Guidance Systems or their respective owners.

Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*

#### Growth Factors

- Conventional (unformulated)
- PODS® - Sustained release

#### Exosomes

- Exo-spin™ - Purification
- ExoLISA™ - ELISA-like detection
- Instant Exosomes™ - purified and characterized
- NTA Service
- Freeze drying service

#### PeptiGel®

- Tunable self-assembling peptide hydrogels

#### Other products and services

- Small Molecules
- Softwell™ - 2D hydrogel (Europe only)
- Orangu™ - Cell counting reagent
- LipoQ™ - Lipid quantification assay
- Primary Hepatocytes

#### Cytogenetics

- Karyotype Analysis
- Array Hybridization

Scan for PODS product page



General [info@cellgs.com](mailto:info@cellgs.com)  
Technical Enquiries [tech@cellgs.com](mailto:tech@cellgs.com)  
Orders and Quotes [order@cellgs.com](mailto:order@cellgs.com)

[www.cellgs.com](http://www.cellgs.com)

**EUROPE**  
**Cell Guidance Systems Ltd**  
Maia Building  
Babraham Bioscience Campus  
Cambridge  
CB22 3AT  
United Kingdom  
T +44 (0) 1223 967316  
F +44 (0) 1223 750186

**USA**  
**Cell Guidance Systems LLC**  
Helix Center  
1100 Corporate Square Drive  
St. Louis  
MO 63132  
USA  
T 760 450 4304  
F 314 485 5424