

# PeptiGel<sup>®</sup> Protocol:

## How to Process PeptiGel<sup>®</sup> Cultures for Cell Staining

---

### Before you begin, please note:

This protocol describes how to process PeptiGel<sup>®</sup> cultures for cell staining including histology, immunohistochemistry and immunofluorescence.

Two options, wax embedding and freezing for cryo-sectioning, are outlined here.

It's recommended to use Greiner Bio-One ThinCerts<sup>™</sup> for this protocol. These inserts are made with polystyrene and will dissolve during the xylene wash step.


This protocol has been written as a guide. Therefore, further assay optimization may need to be carried out for your cultures.

---


### Protocol 1: Wax Embedding

#### Fixation:

1. Prewarm 10% neutral buffered formalin (NBF) solution to room temperature. Immerse the PeptiGel<sup>®</sup> (+/- cells) into the solution and incubate for 30 minutes at room temperature (~25°C).



2. Replace the NBF solution with phosphate buffered saline (PBS) solution, enough to immerse the PeptiGel<sup>®</sup>. Store at 4°C until processing in wax. The maximum recommended storage time is 2 weeks.



#### Processing in Wax:


3. Place the fixed insert containing the PeptiGel<sup>®</sup> cultures (+/- cells) in a cassette. Then place the cassette in a processor basket and assemble into the processor.




- Run the wax processing according to the guideline in the table below.

Reagent	Time (seconds)
50% IMS	90
70% IMS	60
99%IMS	60
99%IMS	60
99%IMS	60
99%IMS	60
99%IMS	60
Xylene	90
Xylene	90
Xylene	90
Molten-wax	90
Molten-wax	120

- Embed the sample according to your instrument’s instructions.





- Cryo-section according to your instrument’s instructions. Individual sections of 5-7 μm thickness are recommended.




## Protocol 2: Freeze for Cryo-Sectioning

### Fast Freeze Using Isopentane:


- Decant liquid nitrogen (LN) into a wide mouthed vacuum flask. Place 25-50 mL of isopentane in a small thick-walled glass beaker and immerse in LN to cool. Lift out when isopentane starts to freeze around the rim of the beaker.
- Apply a drop of optimal cutting temperature (OCT) compound on to a 20 mm cork disc. Place a PeptiGel® (+/- cells) sample onto the disk and apply more OCT until the sample is covered.



3. Place PeptiGel® (+/- cells) with the disc into the isopentane for 5-15 seconds to freeze, then remove.




4. Wrap the frozen PeptiGel® in aluminium foil and store at -80°C.



### Sectioning Using Cryostat:

5. Cryosection according to the instrument's instructions. Cutting sections of 7-10 µm or cutting a 20 µm thickness if culture samples lack structure, is recommended.



For further support, please contact our technical support team at [tech@cellgs.com](mailto:tech@cellgs.com).



General [info@cellgs.com](mailto:info@cellgs.com)  
Technical Enquiries [tech@cellgs.com](mailto:tech@cellgs.com)  
Orders [order@cellgs.com](mailto:order@cellgs.com)

[www.cellgs.com](http://www.cellgs.com)

**EUROPE**  
**Cell Guidance Systems Ltd**  
Maia Building  
Babraham Bioscience Campus  
Cambridge  
CB22 3AT  
United Kingdom  
T +44 (0) 1223 967316  
F +44 (0) 1223 750186

**USA**  
**Cell Guidance Systems LLC**  
Helix Center  
1100 Corporate Square Drive  
St. Louis  
MO 63132  
USA  
T 760 450 4304  
F 314 485 5424