

Matrigen Softwell®

Frequently Asked Questions



Contents

General Questions	3
Cell culture and experimental questions	6
Quality control technical questions.....	8
SoftTrac™ Products.....	9

Matrigen Softwell® Hydrogels

General Questions

- **What is Softwell®?**

Softwell® is a transparent hydrogel bound to the surface of any cell culture well. It is much softer than tissue culture plastic. Softwell® replicates a broad range of tissue stiffness, so you can culture neurons in an environment as soft as the brain, cardiomyocytes on the stiffness of muscle, or whatever your cells' tissue of origin may be. Importantly, this platform can be adapted to most cell-based assays and detection methods, so you can make cell culture more realistic without making it more complicated.

- **How can I order?**

Matrigen products can be ordered in Europe on Cell Guidance Systems website www.cellgs.com or by sending an e-mail to order@cellgs.com. All the products are made to order, which means delivery will typically take 2 – 3 weeks. In other regions, please contact info@matrigen.com.

- **Do you provide smaller pack sizes for initial testing?**

Matrigen products are available for sale as single plates or as packs of 10 dishes/plates.

- **What is the composition of the gels?**

Softwell® gels are made of polyacrylamide crosslinked with bisacrylamide. Polyacrylamide is a non-toxic polymer which forms a soft gel when hydrated. The physical properties of these hydrogels can be modified allowing to control the stiffness. Moreover, polyacrylamide gels may be chemically conjugated with ECM ligands to enable cell attachment, growth and differentiation.

- **Are the gels manufactured under sterile or aseptic conditions?**

All products are manufactured and packaged under aseptic conditions. Plates are individually wrapped in foil pouches.

- **What is the storage and shelf life of the products?**

All products are individually packed and labelled with a manufacturing date. As long as the hydrogels remain sealed or hydrated, all products are stable for 6 months from the date of receipt. Non-Activated and Easy Coat™ recommended storage temperature is 4 - 25°C. Collagen coated hydrogels should be stored at 4°C. Hydration is critical, so keep the hydrogels in the sealed package until ready for use. Adhesion Free products may be stored at room temperature and do not require hydration.

- **Is it possible to open the package and use different wells at different times?**

Yes, this should be possible. To make sure that the hydrogels are kept hydrate, the wells should be filled with PBS until you are ready to use them.

- **What is the range of available stiffness values?**

Softwell® gels are available in 9 standard stiffness values (0.2, 0.5, 1, 2, 4, 8, 12, 25, and 50 kPa), ultra-soft range (0.03, 0.05, and 0.1 kPa) and screen option (range of standard stiffness values).

- **What are the stiffness values for plastic and glass?**

The stiffness values for polystyrene and glass are about 3,000,000 kPa and 60,000,000 kPa, respectively, which makes these materials much stiffer than the hydrogels.

- **Is it possible to have gels with other stiffness values?**

Different stiffness values can be manufactured upon request. An extra charge for custom products will be applied.

- **What formats are Softwell® gels available in?**

Softwell® hydrogels are available in multi-well polystyrene plates, multi-well glass bottom plates, polystyrene dishes, glass bottom dishes, and standard plastic plates with removable coverslips.

- **Which coatings are available in Softwell® gels?**

Softwell® gels are available as Non-Activated, Easy Coat™, Collagen and Adhesion Free™. Non-Activated hydrogels can be used as low protein and cell attachment surface or alternatively as a substrate that can be chemically activated by your own methods. Easy Coat™ hydrogels are activated to bind to any extracellular matrix proteins. Collagen hydrogels are a ready to culture option, since the hydrogels are bound to type I collagen. Adhesion Free™ is a special surface treatment option in which hydrogels resist protein adsorption and are 100% non-adherent to cells.

- **What is the chemistry used on Easy Coat™ hydrogels?**

The exact chemical treatment used on Easy Coat™ hydrogels is proprietary. But essentially, the treatment populates the hydrogel (or surface) with quinone groups, which form covalent bonds with primary amines, thiols, or strong nucleophiles. This means that the activated hydrogels will bind to any peptide or protein. The hydrogels are ready to coat with the ECM protein of choice, e.g. collagen (types I, III and IV), fibronectin, laminin, and poly-lysine. Be aware that once all the binding sites are reacted, additional molecules will not bind directly to the hydrogel.

- **What is the Collagen used for Collagen coated-dishes?**

All collagen-coated products use rat tail collagen type I. Bovine skin type I collagen can be used upon request. All hydrogels are incubated with a solution of 0.02 mg/ml with a constant surface density of collagen across different stiffness values. The exact concentration of collagen in all hydrogels has not been measured.

- **Do you have any experience using gels with other ECM proteins besides Collagen I?**

Based on customers' experience, Softwell® gels bind well to different ECM proteins, including fibronectin, laminin, collagen IV, among others. The gels bind quite readily to the first protein that they are exposed to. Please note that binding normally occurs at a pH range from 6 to 9. We recommend the use of PBS for resuspending the ECM protein of interest.

- **What are the exact dimensions of the plates provided?**

For setting up some experiments, it is important to understand the exact dimensions of the plates provided by Matrigen. If this information is required, please send an e-mail to info@cellgs.com. The standard glass thickness used is #1.5, but #0 can be used upon request.

- **Is it possible to transfer single coverslips to other multi-well plates without damaging the hydrogel?**

The coverslips are smaller than the respective well size and can be carefully removed, while immersed, by using fine-tipped tweezers.

- **What is a High-Throughput Screen (HTS) plate?**

The HTS plates are designed to screen the standard stiffness values. Each row will have the following stiffness values: 0.2, 0.5, 1, 2, 4, 8, 12, 25 and 50 kPa, and three wells with no hydrogel. The wells with no hydrogel will have the same treatment of the hydrogels, which can either be Easy Coat™ or Collagen.

- **Is Softwell® like Matrigel™?**

Hydrogels <1kPa are as soft as Matrigel™. But whereas Matrigel™ is comprised of tumor-derived proteins and growth factors, Softwell® is based on a synthetic polymer.

- **Will the hydrogels detach from the well?**

Hydrogels should remain bound under most conditions, including long-term cell culture. Extra care should be taken during pipetting and aspiration steps for hydrogels <2 kPa.

Cell culture and experimental questions

- **Are the gels safe for my cells?**

Matrigen products are certified non-cytotoxic and pyrogen free.

- **Are the gels compatible with my current assays?**

Softwell® gels are compatible with conventional biological assays and detection methods. The gels are thin allowing easy washing away of reagents and antibodies. Moreover, the gels are optically transparent, meaning that the gels will not interfere with phase contrast or fluorescence microscopy.

- **How do I remove my cells from the gels?**

Cells can be detached and passaged with standard reagents and protocols, just as on tissue culture plastic. The gel is bound to the substrate and will resist detachment.

- **Can I detach my cells with trypsin?**

Yes, but Matrigen recommends incubate the hydrogel in serum-free media or buffer for ≥30 minutes at 37°C to remove serum absorbed in the gel prior to adding trypsin. To bypass this step, use TryPLE™ (ThermoFisher Scientific), which dissociates cells in serum-containing conditions.

- **Following trypsinization, I still cannot extract my cells from the hydrogel. What should I do?**

Some cell types are difficult to extract as they attach strongly to the gel matrix. If the use of TryPLE™ is not sufficient, we recommend the use of the stronger TryPLE™ Select 10X (ThermoFisher Scientific).

- **Can cells be scraped from the gels?**

A cell scraper with silicone rubber blade will do the trick (available from CytoOne®). Scrape with gentle, short strokes to avoid disrupting the hydrogel.

- **Can I fix my cells?**

Softwell® is compatible with fixatives such as formalin, glutaraldehyde, and methanol. Note that dehydrating agents will cause the gel to turn opaque, but regains transparency upon hydration.

- **Is there a recommended protocol for the coating of Easy Coat™ hydrogels with ECM proteins?**

We recommend the use of PBS for resuspending the ECM protein of interest. For coating the hydrogels, we suggest diluting collagen I in PBS at a concentration of 10 µg/ml. For other proteins, the recommended concentration will be much lower (e.g. fibronectin and laminin would be 100-1000 ng/ml). Volume per well would be enough to cover the gel. Finally, a 30 minutes incubation at room temperature or longer is recommended.

- **Will collagen solubilised in acetic acid be detrimental to the gel?**

Acidic solutions will not affect Softwell® gels.

- **What is the recommended product for immunofluorescence staining?**

For immunofluorescence staining, the Softslip™ (hydrogels bound to glass coverslips in multi-well plate) products are ideal since the coverslips can be removed from the well and inverted on a glass slide. Hydrogels bound to glass bottom dishes and 96-well glass bottom plates are also suitable, but keep in mind that the thickness of the gel (60 to 500 µm) will limit the resolution for 40x and 63x objectives. Hydrogels bound to polystyrene multi-well plates and dishes are also compatible for immunofluorescence, but the thickness of the plastic and gel will also limit the resolution for 40x and 63x objectives.

- **Can I image my cells?**

As mentioned above, the hydrogels are transparent and compatible with brightfield, phase contrast, and fluorescence microscopy. The combined thickness of the hydrogel and substrate is within the focal range of 40X objectives. For higher magnifications, use Softslip™ hydrogels, which can be inverted onto a coverslip.

- **Should I expect auto-fluorescence from Softwell® gels or from antibodies binding to the gels when performing immunofluorescence assays?**

Softwell® gels do not have any associated auto-fluorescence and they will not bind to anything besides the ECM proteins used for cell adherence. To ensure removal of unbound detection molecules, perform wash steps 3x longer than directed by standard protocols. It is also recommended to use 1% BSA and/or 5% goat serum as a blocking buffer to minimize any background in the images.

- **Do Softwell® gels keep a consistent focal plane?**

Softwell® hydrogels have a uniform surface, so a constant focal plane can be maintained if the plate is represented by a single stiffness. However, more adjustments might be required with variable stiffness plates, because the hydrogels swell to a different extent according with their stiffness (softer gels swell more than stiffer gels).

- **How do I isolate RNA?**

Matrigen recommends using a spin-based kit such as RNeasy® (Qiagen). Incubate the hydrogel with lysis buffer for 10 minutes on ice, collect the lysate, and repeat to recover RNA absorbed in the gel.

Quality control technical questions

- **What is the methodology for stiffness calibration per batch?**

Stiffness calibration and Young's Modulus calculus are performed by applying Hertz's sphere model (assuming a Poisson ratio of 0.48 for polyacrylamide) after the measurement of indentation of a < 1 mm carbide steel sphere (density = 15.63 g/cm³) upon the surface of the hydrogel. The measurements are taken while the hydrogel and sphere are immersed in phosphate-buffered (PBS) (density = 1 g/cm³). This procedure is applied to each batch, meaning that each time a stock (prepolymerization) solution of a targeted stiffness is made, an aliquot is polymerized for stiffness validation by using the above described method.

- **Is stiffness calibration performed with Atomic Force Microscopy (AFM)?**

Stiffness calibration is performed by applying Hertz's theory as described in the previous question. The obtained values are well correlated with AFM measurements. AFM thickness validation is not routinely performed.

- **What is the calibration for gel thickness?**

The gel thickness may vary depending on stiffness and format within a range of 200 and 500 µm. For hydrogels on polystyrene, the thickness is ~ 500 µm, whereas for hydrogels on Softslip™, 6, 12, 24, and 96 well glass bottom plates, the thickness is ~200 µm. Finally, for hydrogels on Softview™, the thickness will be ~100-150 µm. The gels are no thicker than 500 µm to minimize the interference of the gel in microscopy and other assays. The gel thickness can be modified upon request. The stiffness of the hydrogel does not change after coating.

- **What is the size of the pores in the gels?**

All the pores in the gels are <0.2 µm. To ensure that this parameter is satisfied, the gels are overlaid with 0.2 µm fluorescent beads to confirm that they cannot penetrate the softest gels.

- **Do you perform electron microscopy to analyze gels morphology?**

Matrigen has not performed electron microscopy analysis on the gels.

- **What is the methodology for quality and uniformity control of the collagen coated gels?**

The surface density of collagen is relatively constant across stiffness. This uniformity was confirmed by using chemiluminescence-based detection of gel-bound collagen across a range of different stiffness values.

SoftTrac™ Products

- **What are SoftTrac hydrogels?**

SoftTrac™ are hydrogels with fluorescent microspheres immobilized at the surface.

- **What will be the use of SoftTrac hydrogels?**

SoftTrac™ hydrogels main use is for traction force microscopy.

- **What microspheres are immobilized in the hydrogels?**

Matrigen provides hydrogels with the following range of microspheres: 0.2 µm yellow/green, 1 µm yellow/green, and 1 µm red microspheres. Yellow-green fluorescent microspheres (0.2 or 1 µm diameter) with excitation and emission wavelengths of 505 and 515, respectively. Red fluorescent microspheres (1 µm diameter) with excitation and emission wavelengths of 580 and 605, respectively. The density of the beads can be adjusted upon request.

Cell Guidance Systems' reagents and services enable control, manipulation and monitoring of the cell, both *in vitro* and *in vivo*.

Growth Factors

- Recombinant
- Sustained Release

Exosomes

- Purification
- Detection
- Tracking
- NTA Service

Small Molecules

Cell Counting Reagent

Matrix Proteins

Cell Culture Media

- Pluripotent Stem Cells
- Photostable
- *In Vitro* Blastocyst Culture
- ETS-embryo Culture
- Custom Manufacturing Service

Gene Knock-Up System

Cytogenetics Analysis



General info@cellgs.com
Technical Enquiries tech@cellgs.com
Quotes quotes@cellgs.com
Orders order@cellgs.com

www.cellgs.com

EUROPE
Cell Guidance Systems Ltd
Maia Building
Babraham Bioscience Campus
Cambridge
CB22 3AT
United Kingdom
T +44 (0) 1223 967316
F +44 (0) 1223 750186

USA
Cell Guidance Systems LLC
Helix Center
1100 Corporate Square Drive
St. Louis
MO 63132
USA
T 760 450 4304
F 314 485 5424